

URIT

Compact design and easy operation

URIT-50 Urine Analyzer



Global Diagnostics Supplier

URIT-50

Urine Analyzer

FDA 510 (K) No.K082811 Approved



Specifications

Parameters: pH, Nitrite, Specific Gravity, Ascorbic Acid, Blood, Glucose, Bilirubin, Urobilinogen, Ketone, Leucocytes, Protein
Measurement Principle: Dual wavelength reflectance method
Measuring Wavelength: 3 LED Wavelength (550,620,720nm)
Sample supply: manual dipping
Reaction time: 60 seconds
Throughput: Normal mode, 60 samples /hour Continuous mode, 125 samples /hour
Display: LCD screen
Memory capacity: 5000 results
Strip: URIT 10G /11G urine reagent strips
Interface: Serial port(RS232), parallel port (For external printer), PS/2(For barcode reader)
Work condition: Environmental Temperature: 15-30°C, RH: 20%~80%
Power supply: Adopter 12V DC, 4.0A
Dimension: 290mm x 200mm x 100mm
Weight: 2kg



URIT URINE CONTROL

Features

- ✔ Cold light source facility with long life and high stability
- ✔ Internal thermal printer
- ✔ Easy operation, cost-efficiency and high quality
- ✔ Rs232 serial interface to host PC/LIS, Barcode reader optional
- ✔ URIT UQ-10 & UQ-11 Urinalysis Control



URIT 10G、11G urine reagent strip

Complete urine analysis with Ascorbic Acid

URIT MEDICAL ELECTRONIC (GROUP) CO., LTD.

Head Office

Add.: No.4 East Alley, Jiuhua Road, Guilin,
Guangxi 541001, P.R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

E-mail: export@uritest.com service@uritest.com

Shenzhen Office

Add.: Room EJ, 24th Floor, New Baohui Building, Nanhai Road,
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150

Http://www.urit.com



ISO9001
ISO13485